









Table 1: Tabulated MTI Mechanical Convection Oven

Model	Chamber Dimension L x W x H	Volume	Working Temperature	Features	Power	Item Image
DHG-9015	10" x 10" x 8"	11.5L	Max: 300°C (<30 mins) Cont: 250°C	Convection Drying Oven with Digital Temperature Controller	800W	
BPG-7032	13" x 12" x 12"	30L	Max: 250°C	NRTL/CSA Certified Convection Oven with Digital Temperature Controller	850W	
BPG-7082	17.7" x 15.7" x 17.7"	80L	Max: 250°C (slow heated) Cont: 200°C	NRTL/CSA Certified Convection Oven with Digital Temperature Controller	1500W	

EQ-DHG-9140A	21.6" x 21.6" x 17.7"	136L	Max: 250°C	Convection Drying Oven with Digital Temperature Controller	2300W	
EQ-DHG-9245A	23.5" x 17" x 29.5"	190L	Max: 300°C (< 30 mins) Cont: 250°C	Convection Drying Oven with Digital Temperature Controller	2450W	
EQ-DHG-9425	25" x 23" x 53"	420L	Max: 250°C	Large Mechanical Convection Oven with Digital Temperature Controller	3500W	
EQ-DHG-9000JB	14" x 14" x 14"	43L	Max: 400°C	CE Certified Mechanical Convection Oven with 28 Segments Programmable Controller	2200W	

KSL-500X-71	335mm x 450mm x 450mm	67L	Max: 500C	Convection oven with Programmable Temperature Controller	2200W	
-------------	-----------------------	-----	-----------	--	-------	---