






Table 1: Tabulated MTI Vacuum Oven

Model	Chamber dimension L x W x H	Volume	Working temperature	Features	Power	Item image
DZF-6020- HT400P	12" x 12" x 11"	25L	Max: 400°C	Vacuum Oven with Quartz Window, Flow-meter & Water Chiller	1600W	
DZF-6050-HT	16" x 13" x 14"	53L	Max: 500°C (< 8 hours) Cont: 450°C	30 Segments Temperature Controller, Flowmeter, Chiller and Vacuum Pump System	3500W	
DZF-6090-HT	18" x 18" x 18"	91L	Max: 550°C (< 2 hours) Cont: 500°C	Floor-Stand Vacuum Oven with 30 Segments Temperature, Two Heating Zones and Vacuum Pump & Water Chiller	2600W	
DZF-6090-II	18" x 18" x 18"	91L	Max: 200°C	Vacuum Oven with Vacuum Pump and Two Heating Chambers	2200W	

DZF-6210	22" x 25" x 24"	215L	Max: 200°C	Large Floor Stand Vacuum Oven with Vacuum Pump and Three Heating Zone	2050W	
DZF-3120	800mm x 500mm x 300mm per level	360L	Max: 150°C	Larger Capacity Vacuum Drying Oven with Tri-level Shelf Heating Modules	18000W	